

# Foot Machine



User Manual

# Foot Massager

Thank you for purchasing For proper operation and use of this foot massager, please read the operating manual carefully and pay special attention to safety.

This manual provides safe and effective operation. Please do not lose it for future reference.

---

## CATALOG

Safety guidelines . . . . . 2 - 5

Product operation and product maintenance.....6

Troubleshooting Guide..... 7

Warranty and Specification Information.....8

# Safety guidelines

---

## ( Please consult the doctor before using )

- 1a. Those who are undergoing medical treatment or surgery are not recovered.
- 1b. Patients with hypertension, heart disease, cerebral arteriosclerosis and cerebral thrombosis.
- 1c. Patients with high fever, malignant tumor and osteoporosis  
Pregnant women during menstrual period
- 1d. A cardiac pacemaker or foot implant device is being used.  
Other acute patients.
2. This product has surface heating function and should be noted for heat sensitive people.
3. Children should be taken care of to avoid playing with toys.
4. Children aged 12 or above can use this product. When using this product, disabled persons or persons with restricted behavior should be guided and used by caregivers.
5. This product can not be used for children or people in a state of unconsciousness. And need regular inspection and addition.  
Thermal function contacts the body position to avoid scalding.
6. The foot machine must have a socket fitted with the plug.  
Improper connection of the plug may cause electric shock. If you have any doubts about the safety of this product or the socket is not suitable for the plug of this product, please contact the professional and technical personnel for inspection or installation.

## Safety guidelines

---

7. Products must be kept away from sharp objects, such as damaged cloth. In order to avoid toes involved in the machine cause personal injury.
8. If the product falls into the water accidentally, it must be maintained and not allowed to be used.

### **Use of the environment**

1. Do not use the product in wet environment such as bathroom or near water source in order to avoid malfunction and even endanger personal safety.
2. Do not use this product at room temperature higher than 40 degrees Celsius or lower than 5 degrees Celsius, because the thickening of lubricating oil may increase the wear of machine gear transmission parts and affect the service life.
3. Please do not use in the open air or for prolonged exposure to the hot sun, which may damage the paint layer.
4. If the product is transferred to warm environment at low temperature, the internal metal parts may gather water vapor and affect the normal operation of the product, even cause damage to the machine. It is suggested that the product be placed at normal room temperature for one hour before starting to use after transferring the environment.

# Safety guidelines

---

## **Precautions**

1. Do not pull the power cord or use the power cord to suspend the product, so as not to cause wire breakage or bad contact. When pulling out the power plug, you should hold the plug and pull it out.
2. If the product power cord is damaged, it must be replaced by a professional
3. Do not cover the blanket on the product or block the heat sink under the machine when using, so as not to cause the temperature of the internal parts to be too high and lead to failure.
4. When using, please press the switch lightly, avoid force or nail operation.
5. Do not allow the control panel to contact liquids such as water.
6. Do not insert or drop objects into any openings of the product. Do not sleep while using products.

# Safety guidelines

---

## Precautions

7. Do not use this product for one hour before and after eating, or for more than 30 minutes continuously, in order to avoid adverse effects. Do not pull the power cord or use the power cord to suspend the product, so as not to cause wire breakage or bad contact. When pulling out the power plug, you should hold the plug and pull it out.  
from the socket before use or cleaning.
8. Do not use this product in drunkenness or untimely use.
9. If you feel unwell during use, you should stop using it immediately.
10. If you want to cut off the power supply, please close all the operation control and pull the power plug from the socket.
11. Please remove the power plug from the socket before use or cleaning.
12. If the machine fails, it should be stopped immediately and feed back the problem to our staff
13. This product can not be used for self-diagnosis and treatment. If it needs diagnosis or treatment, it should be handled by professional medical personnel.

## Product operation

---

Power switch: press the power switch and the machine starts to work.



Scraping: two kinds of foot scraping mode selection, but also through this button to scrape foot function closed.



Air pressure: two air pressure intensity selector, and the pressure function can also be closed by this button.



Warm feeling: press the heating button, the heating function is turned on, or the heating function can be closed by this button.



Negative Ion: Press the negative ion button, and the negative ion function can be turned on, or it can be turned off by this button.

## Product maintenance

---

1. Storage : When not in use, put the product in a color box and place it in a cool place. Do not place products in direct sunlight or high temperature.
2. Cleanliness : Before cleaning, confirm that the main power switch is closed and pull out the power plug. Please use dry cloth to clean the surface of the product. Do not use corrosive cleaners, such as benzene or thinner, to clean products.
3. Maintenance : If you do not use this product for a long time, please use a dust cover to prevent dust falling.

# TROUBLESHOOTING GUIDE

---

If problem are not from above reason, please turn off power supply and contact the distributor or manufacturer for assistance. Note: Do not dismantle the produce by yourself, and our company shall shoulder no responsibility to the damage of the product or the inquiry of people caused by improper use. Thanks for your cooperation!

Problem / Issue	Possible Cause	Method
No operation	<ol style="list-style-type: none"><li>1. Non electricity</li><li>2. The function key is not opened.</li></ol>	<ol style="list-style-type: none"><li>1. Access to power or check whether the power is connected properly.</li><li>2. Turn on on / off button.</li></ol>
Sudden stop operation	<ol style="list-style-type: none"><li>1. Power off</li><li>2. Massage time is more than 15 minutes</li><li>3. Continuous use time is too long, product overheating protection</li><li>4. Excessive force, massage ball blockage</li></ol>	<ol style="list-style-type: none"><li>1. Check whether the power interface is properly connected.</li><li>2. Reboot</li><li>3. Restart after cooling</li><li>4. Reboot</li></ol>



# WARRANTY

---

1. Warranty rights are not transferable.
2. The product warranty is valid for one year from the date of purchase.
3. You can enjoy the free warranty during the warranty period (the appearance of the product is not covered by the warranty).
4. The product warranty requires proof of purchase and warranty card.
5. The above warranty does not apply to commercial users such as leasing, trading, and institutions.
6. The following conditions are not covered by the warranty, only promised to provide paid repair services (changes for labor and materials).
  - a. Damage caused by all human factors: including use in abnormal working environment, use and storage in accordance with the instructions, damage caused by improper use, etc.;
  - b. The user or the third party not authorized by the company to dismantle, modify and repair.
  - c. Damage caused by transportation after purchase.
  - d. Damage caused by other force majeure factors such as floods, fires, lightning strikes and earthquakes.

## SPECIFICATION INFORMATION

Product Name		Two Face Cavalier	
Model no.	OS-808	Colour	Black / White
Rated Voltage	100V~240V	Rated Frequency	50Hz
Rated Power	40W		